



Live Freely. Breathe Confidently.

The Kingon P2-E gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators.

Whether relaxing at home, heading out for daily activities, or travelling, the P2-E provides reliable oxygen therapy in a highly compact form.

Simple, User-Friendly Design

The Kingon P2-E features an easy-to-navigate interface with a bright 2.8" LCD screen, ideal for users with impaired vision. Adjustable flow settings allow the device to adapt as your oxygen needs change. Visual alerts remain on-screen even when silenced, ensuring important information, like low battery, low oxygen purity, or no breath detected - is always visible.



Reliable Performance

Ideal for home and on-the-go oxygen therapy.



Simple Maintenance

Easy-to-clean filters and accessible parts.



Quiet Operation

Only 37 dB, as quiet as library ambience.



Trusted Support

Comprehensive Australian warranty and service.

P2-E Portable Oxygen Concentrator

Lightweight power. Reliable oxygen.
Everyday freedom.

The Kingon P2-E is the world's smallest and lightest portable oxygen concentrator in its class. Weighing just 1.98 kg and measuring 22.1 × 8.5 × 16 cm. Compact and easy to carry, it can be lifted effortlessly with one hand, whether you are at home or on the go.



-  **Lightweight 1.98 kg** - compact and portable
-  **Up to 5 hours battery life** - stay powered longer
-  **Smart pulse-dose delivery** - oxygen only when you inhale
-  **Quiet 37 dB operation** - low noise for discreet use
-  **Airline approved** - travel confidently

At a Glance

Specification	Details
Oxygen Concentration	90% (−3% / +6%) at all settings
Flow Settings	1–5 (pulse dose delivery)
Weight	1.98 kg (with battery)
Size	22.1 L × 16 H × 8.5 W cm
Noise Level	37 dB (setting 2)
Battery Duration	1: 5h, 2: 4h, 3: 3h15, 4: 2h30, 5: 1h42
Recharge Time	Approx. 4 hours (AC/DC)
Power Options	AC 100–240 V · DC 12 V · Battery
Display	2.8" LCD screen
Warranty	Device: 3 yrs standard + 2 yrs extension upon registration · Accessories: 12 mths

Package Inclusions

- | | | | |
|---|---------------------------------------|---|-------------------|
|  | 1 x Kingon |  | 1 x Carry Bag |
|  | 2 x Batteries |  | 1 x Nasal Cannula |
|  | 1 x AC Charger |  | 5 x Air Filters |
|  | 1 x DC Charger
(charge in the car) | | |

Discover the Kingon difference.

Visit cpapclub.com.au to learn more.