

P2-E6 Portable Oxygen Concentrator

Lightweight power. Reliable oxygen.
Everyday freedom.

The Kingon P2-E6 is the world's most versatile portable oxygen concentrator, combining lightweight design with high-performance oxygen delivery. Weighing just 1.98 kg and measuring 22.1 x 8.5 x 16 cm, it's engineered for mobility, whether at home or on the go. Lift it with ease or carry it in its compact bag for added convenience.



-  **Lightweight 1.98 kg** - compact and portable
-  **5h 24min battery life** - stay powered longer
-  **Smart pulse-dose delivery** - oxygen only when you inhale
-  **Quiet 37 dB operation** - low noise for discreet use
-  **Airline approved** - travel with confidence

At a glance

Specification	Details
Oxygen Concentration	90% (-3% / +6%) at all settings
Flow Settings	1 - 6 (pulse dose delivery)
Weight	1.98 kg (with battery)
Size	22.1 L x 16 H x 8.5 W cm
Noise Level	37 dB (setting 2)
Battery Duration	1: 5h 24m 2: 4h 30m 3: 3h 18m 4: 2h 30m 5: 2h 06m 6: 1h 42m
Recharge Time Power	Approx. 4 hours (AC/DC) AC 100 - 240 V · DC 12 V Battery
Options Display	2.8" LCD screen
Warranty	(3-year standard + 2-year free extension when you register your device) 1 year for battery & sieve bed

P2-E6 Package Includes

-  1 x Kingon P2-E6
-  1 x Carry Bag
-  2 x Batteries
-  1 x Nasal Cannula
-  1 x AC Charger
-  5 x Air Filters
-  1 x DC Charger (charge in the car)

Discover the Kingon difference.

Talk to your oxygen provider or visit kingon.com.au to learn more.



Live Freely . Breathe Confidently .

The Kingon P2-E6 gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators.

Whether relaxing at home, heading out for daily activities, or travelling, the P2-E6 provides reliable oxygen therapy in a compact, easy-to-use form.

Simple, User-Friendly Design

Featuring an easy-to-navigate interface with a bright 2.8" LCD screen, the P2-E6 is ideal for users with impaired vision. Adjustable flow settings allow the device to adapt as your oxygen needs change. Visual alerts remain on-screen even when silenced, ensuring important information - like low battery, low oxygen purity, or no breath detected - is always visible.



Reliable Performance

Perfect for home and on-the-go supplemental oxygen therapy



Simple Maintenance

Easy-to-clean filters and accessible parts.



Quiet Operation

Only 37 dB - as quiet as a soft conversation.



Trusted Support

Comprehensive Australian warranty and service.

Discover the Kingon difference.
Visit cpapclub.com.au to learn more.