

P2-E7 Portable Oxygen Concentrator

Lightweight power. Reliable oxygen.
Everyday freedom.

The Kingon P2-E7 is the world's smallest and lightest portable oxygen concentrator. Four Cylinder Compressor. Weighing just 1.98 kg and measuring 22.1 × 8.5 × 16 cm, it's the closest you can get to a pocket-sized POC. Take it anywhere—at home or on the go—and lift it easily with one



- Lightweight 1.98 kg** - compact and portable
- Up to 7 hours battery life** - stay powered longer
- Smart pulse-dose delivery** - oxygen only when you inhale
- Quiet 37 dB operation** - low noise for discreet use.
- Airline approved** - travel confidently

At a Glance

Specification	Details
Oxygen Concentration	90% (−3% / +6%) at all settings
Flow Settings	1–7 (pulse dose delivery)
Weight	1.98 kg (with battery)
Size	22.1 L × 16 H × 8.5 W cm
Noise Level	37 dB (setting 2)
Battery Duration	Up to 7 hours duration on setting 1
Recharge Time	Approx. 4 hours (AC/DC)
Power Options	AC 100–240 V · DC 12 V · Battery
Display	2.8" colour LCD screen
Warranty	3 yrs device + 2 yrs extended on registration of device. · 12 mths sieve bed/ batteries

Package Inclusions



Discover the Kingon difference.

Talk to your oxygen provider or visit kingon.com.au to learn more.



Live Freely. Breathe Confidently.

The Kingon P2-E7 gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators.

The Kingon P2-E7 gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators. Whether relaxing at home, heading out for daily activities, or travelling, the P2-E7 provides reliable oxygen therapy in a highly compact form.

Simple, User-Friendly Design

The Kingon P2-E7 features an easy-to-navigate interface with a bright 2.8" colour LCD screen, ideal for users with impaired vision. Adjustable flow settings allow the device to adapt as your oxygen needs change. Visual alerts remain on-screen even when silenced, ensuring important information—like low battery, low oxygen purity, or no breath detected—is always visible.



Reliable Performance

Ideal for home and on-the-go oxygen therapy.



Simple Maintenance

Easy-to-clean filters and accessible parts.



Quiet Operation

Only 37 dB — as quiet as library ambience.



Trusted Support

Comprehensive Australian warranty and service.